

Clearing and Staining Protocol - Draft 5

Harold Hodgins

June 15, 2016

The following procedure is a combination of the protocols of Dr. Armbruster [1] and Dr. Helland[2]. It has been modified specifically for use with zebra fish in small petri dishes. The volumes for the stock chemicals and required solutions are based on 4 fish and 10mL / required solution / fish. Rounded to 50mL total volume of each required solution.

1 Stock Chemicals

- Glacial Acetic Acid - 15mL
- Alcian Blue - cheapest source is fine - 10mg
- Alizarin red - a little goes a long way - 0.5mg
- 95% Ethanol - 165mL
- Potassium Hydroxide Crystals - 0.5g¹
- Sodium Borate - borax from grocery store is fine - 80mL saturated solution
- Trypsin Powder - cheapest source is fine - $\approx 0.5g$ total

If preserving

- 10% Formalin - 50mL

If Bleaching

- Hydrogen Peroxide - pharmacy grade (4%) is sufficient, replace every 6 months - 7.5mL

If Storing in Glycerin

- Glycerin - cheapest source is fine
- Thymol Crystals - a little goes a long way

¹If not bleaching. If also bleaching then 1g total. Glycerin storage amounts not included since we aren't storing our specimens in Glycerin

If Storing in Resin

- Clear polyester casting resin
- Catalyst

2 Required Solutions

- Preservation/Dehydration
 - 50mL 10% Formalin
 - 50mL 30% Ethanol
 - * 16mL 95% Ethanol
 - * 34mL DI water
 - 50mL 50% Ethanol
 - * 26mL 95% Ethanol
 - * 24mL DI water
 - 50mL 70% Ethanol
 - * 37mL 95% Ethanol
 - * 13mL DI water
 - 50mL 95% Ethanol
- 50mL Cartilage Stain
 - 15mL acetic acid
 - 35mL 95% Ethanol
 - 10mg Alcian Blue stain
- 50mL Bone Stain
 - 50ml 1% pottassium hydroxide
 - 0.5mg Alizarin Red stain
- Neutralization
 - 50mL saturated sodium borate
- Bleaching
 - 7.5mL 4% Hydrogen Peroxide
 - 42.5mL 1% Potassium Hydroxide
- Clearing Stage 1,2
 - 100mL 30% sodium borate
 - 0.5g trypsin

3 Procedure

1. Obtain the fish post-mortem.
2. Preserve the fish.
 - (a) place in **10% formalin** for 5 days
 - (b) rinse with water
 - (c) store in water for 2 days
 - (d) rinse with water
 - (e) store in **30% ethanol** for 2-5 days (<150 mm 2 days)
 - (f) store in **70% ethanol** for 2-5 days (<150 mm 2 days)
3. **Remove the guts (optional)**
4. Dehydration
 - (a) rinse with fresh DI water
 - (b) transfer to **50% ethanol**
 - 10-20mm : 1 day
 - 20-80mm : 2 days
 - 80-200mm : 3 days
 - >200mm : 5 days
 - (c) place in **95% ethanol**
 - 10-20mm : 1 day
 - 20-80mm : 2 days
 - 80-200mm : 3 days
 - >200mm : 7 days
5. Cartilage Staining
 - (a) place specimens in **30/70 acetic acid/absolute ethanol+ 20mg alcian blue stain (per 100ml)**
 - (<80mm) 1 day
 - (80-500mm) 1.5 days
 - (>500mm) 2 days
6. Neutralization
 - (a) place specimens in **saturated sodium borate** directly post staining
 - (<100mm) $\frac{1}{2}$ day
 - (>100mm) 2 days (may need to change solution half way through)
7. **Bleaching (optional)**

- (a) place specimens in **15/85 4% hydrogen peroxide/1% potassium hydroxide**, most will depigment within an hour, darker will take longer, ***DO NOT LEAVE IN TOO LONG***
 - (b) If bubbles form in skeleton immediately move to 1% potassium hydroxide solution
8. Clearing (Stage 1)
- (a) place specimens in **30% sodium borate + trypson powder** for 7-10 days till 60% clear (you can see the vertebral column but no fully clear yet)
 - (b) change solution every 10 days
 - (c) change solution when it turns blue
 - (d) keeping specimens in light may help speed things up
9. Bone Staining
- (a) place specimens in **1% potassium hydroxide + 1mg Alizarin Red**
 - (b) leave till bones are desired darkness level
 - 10-80mm : 1 day
 - 80-200mm : 2 days (replace solution after 1 day)
 - > 200mm : 4 days
 - (c) ***DO NOT LEAVE IN TOO LONG***
 - (d) Helland recommends staining with a stirrer at room temp in the dark
10. Clearing (Stage 2)
- (a) place specimens back into **30% sodium borate + trypson powder** for 7-10 days till fully clear
 - (b) change solution when it turns blue
 - (c) keeping specimens in light may help speed things up
11. Preservation
- If in Glycerin
 - (a) put specimens in **30/70 of glycerin/1% potassium hydroxide** for 2-7 days (depending on size)
 - (b) put specimens in solution **60/40 of glycerin/1% potassium hydroxide** for 2-7 days (depending on size)
 - (c) put specimens in **glycerin with thymol** for final storage
 - (d) move to next solution after the specimen(s) sink to the bottom
 - If in Resin (per specimen)
 - (a) Mix correct volume of resin and catalyst
 - (b) Pour in tube around specimen
 - (c) top off resin and let fully cure

4 Time Line

Stage	Time (days)
Aquisition	1
Preservation	11-17
Gutting	1
Deyhydration	2-7
Cartilage Staining	1-2
Neutralization	1-2
Bleaching	1
Clearing 1	7-10
Bone Staining	1-3
Clearing 2	7-10
Preservation	4-14
Total (days)	37-74
Total (weeks)	5-10

References

- [1] Jonathan W. Armbruster. Clearing and staining methods. http://www.cypriniformes.org/files/Clearing_and_staining.pdf.
- [2] Synnove Helland. Diagnostics - staining protocol of cartilage & bone. "www.feap.info/shortcut.asp?FILE=1082".